

PHILIPS
RESPIRONICS

DreamWear

Rediscover dreams

Like wearing nothing at all*

top of head connection

gives your patient the freedom to move

Soft flexible frame

Provides a soft touch on the patient's face

Small	1116745
Medium	1116746
Large	1116747



Fabric wraps

Provides a soft interface to the patient's cheek

Wraps (x2)	1116754
------------	---------

Slim-line headgear

Headgear	1116750
----------	---------

Under-the-nose nasal cushion

Prevents red marks, discomfort or irritation in the nostrils or on the nose bridge**

Under-the-nose cushion

Small	1116740
Medium	1116741
Large	1116742
Medium wide	1116743

Fit Pack includes all 4 cushion sizes (S, M, L, MW) and 1 medium sized frame.

Mask with headgear	International ⁽²⁾	Global ⁽¹⁾
Medium	1116720	1116700

*2015 Philips User Preference Questionnaire, data on file.

**Mask does not directly contact the bridge of the nose or nostrils.

Helpful tips

Here are some tips to keep in mind when fitting and using your DreamWear sleep therapy mask.

To learn more please visit our video guide at philips.com/dreamfamily



- Before using the mask, hand wash the entire mask before use (see Mask Care in your Instructions for use). Wash your face. Do not use moisturizer or lotion on your hands or face.
- Four cushion sizes (S, M, L, MW) are available. Place the cushion-sizing gauge under your nose to determine your size. For the best performance, use the smallest cushion that fits your nose.
- The most common mistake when fitting is over-tightening the headgear. Signs of overtightening include redness, sores, or bulging skin around the edges of the mask.

Fitting instructions



1 Hold the sizing gauge horizontally under your nose and select the cushion size based on where the outer edge of your nostrils and tip of your nose contact the gauge.



2 The cushion opening will be directly under your nostrils. The cushion will hug your nose and leak will be minimal.

(At no time should your nose be inside the cushion opening. If any part of the cushion overlaps your nose, the cushion is being worn incorrectly.)



3 With the mask assembled, place the cushion under your nose.



4 Position the frame on the top of your head.



5 Pull the headgear over the back of your head.



6 To adjust the mask, peel the headgear tabs away from the fabric. Adjust the straps length. Press the tabs back against the fabric to reattach.

(Note: Do not overtighten the headgear.)



7 Position the mask until it fits comfortably. When finished, the elbow should rest at the top of the head.

Mask frame sizing

The medium (MED) mask frame will comfortably fit most faces. If the MED frame does not fit your face, contact your healthcare provider to see if a small (SM) or large (LG) mask frame would better suit your needs.

(1) Global: Brazilian Portuguese, Dutch, English, French, German, Italian, Spanish.

(2) International: Arabic, Bulgarian, Croatian, Czech, Danish, Dutch, English, Estonian, Finnish, French, German, Greek, Hebrew, Hungarian, Italian, Kazakh, Latvian, Lithuanian, Macedonian, Norwegian, Polish, Portuguese, Romanian, Russian, Serbian, Slovak, Slovenian, Spanish, Swedish, Turkish, Ukrainian.

© 2015 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.com/respironics
Customer Service +1 724 387 4000
PN 1126891 MCI 4107166

