



Your Healthy choice  
in CPAP accessory disinfection  
with simple one-touch operation.

## What is Lumin™?

Lumin is the easiest and fastest way to disinfect a CPAP mask and accessories. Lumin works with a disinfection cycle time of 5 minutes, no harmful ozone, and a 99.9% kill rate for harmful bacteria, viruses, mold and fungus. It is the ideal high margin retail accessory for a DME servicing CPAP patients.

### Lumin™

ITEM NO. LM3000  
GTIN: 817783020404  
MSRP \$299.00  
MAP \$249.00





## CPAP MASK WIPES

For a quick and easy pre-treatment, effectively remove dirt, grease, oils and other organic residue from your mask with ease and comfort.



**Aloe**

ITEM NO. M200

GTIN:817783020084



**Citrus**

ITEM NO. M201

GTIN: 817783020091

**3B Medical, Inc.**

Add: 799 Overlook Drive, Winter Haven FL 33884

Tel: 863-226-6285 Fax: 863-226-6284 [www.3Bproducts.com](http://www.3Bproducts.com)

[Info@3BLumin.com](mailto:Info@3BLumin.com)

[www.3BLumin.com](http://www.3BLumin.com)